

## Honey and Biotech

### Is raw honey sterile?

**General:** Honey is the oldest wound dressing material known to humans. Per the National Institute of Health, studies reported that honey has antioxidant, antibacterial and anti-inflammatory properties. It can be used as a wound dressing to promote rapid and improved healing. These effects are due to honey's antibacterial action. Healing is also promoted due to its high acidity, osmotic effect, antioxidant content and hydrogen peroxide content. This experiment initiates the investigation of raw honey and its sterility by using serial dilution and spread plating methodology to enumerate the number of bacteria within a raw sample.

#### Vocabulary

Define the following words in your own terms and cite sources used:

serial dilution

CFU

antibacterial

antioxidant

anti-inflammatory

raw honey

**Problem:** Is raw honey purchased from local markets sterile?

**Hypothesis:**

**Variables:**

**Independent:**

**Dependent:**

**Controls:**

**Safety:** eye protection

#### Materials

nutrient agar plates as prepared by [SLOP#99](#) (in Soy Fresh Unit found at [grownextgen.org](http://grownextgen.org)) or already

prepared plates

distilled water

micro tubes

H1000 micropipette

H200 micropipette

micropipette tips

honey samples

**Tools & equipment:** vortexer, electronic balance, incubator

## Honey and Biotech

### PROCEDURE

#### Prep Honey Sample:


1. Label a micro tube as  $10^0$ ; add initials and honey sample number.
2. Using a H1000 micropipette, add 900 $\mu$ L of distilled water into the labeled micro tube.
3. Using a H200 micropipette, add 100 $\mu$ L of honey sample to the micro tube with distilled water and mix for 15 seconds on a vortexer. (If honey is too thick, warm in a 60 $^{\circ}$ C water bath for 2-3 minutes before pipetting.)
4. Repeat steps 1-3 for each honey sample being investigated.

#### Set Up Serial Dilutions

1. Label micro tubes  $10^1$  through  $10^3$ ; add initials and honey sample number.
2. Use a H1000 micropipette, to add 900 $\mu$ L of distilled water into each micro tube.
3. Then use a H200 micropipette to add 100 $\mu$ L from the  $10^0$  tube to tube labeled  $10^1$  and mix on the vortexer for 15 seconds.
4. Using aseptic technique, transfer 100 $\mu$ L of sample in  $10^1$  tube to the labeled  $10^2$  tube and mix on the vortexer for 15 seconds.
5. Using aseptic technique, transfer 100 $\mu$ L of sample in  $10^2$  tube to the labeled  $10^3$  and mix on the vortexer for 15 seconds.

#### Set Up Spread Plates

1. Along the outer rim of the bottom of each nutrient agar plate, label the plate with corresponding dilution factor, honey number, date and initials.
2. Using aseptic technique, pipet 100 $\mu$ L of each diluted sample in the middle of the corresponding plate.
3. Using a plate spreader or sterile swab, move the sample over the entire plate. If using a sterile swab, use a new one on each plate. If using a plate spreader, be sure to rinse the plate spreader in between each sample.
4. Example diagram for setting up serial dilutions and plating: (adapted from Bio-Rad)



## Honey and Biotech

### Data

Create a data table in APA Format

### Conclusion

REE: (Restate evidence; discuss if hypothesis is accepted, rejected or inconclusive and give actual data evidence to support)

PE: (Possible/Potential errors that occurred during the investigation -- reminder your lab partner is not an error)

PA: (Practical applications; importance of investigation)

## Honey and Biotech

\*This document may be reproduced for educational purposes, but it may not be reposted or distributed without crediting GrowNextGen and The Ohio Soybean Council and soybean checkoff.